



# Support Growth

Fibra-Cel® Disks – A solid support growth material for cell culture





# Suspend Your Disbelief

By using Fibra-Cel Disks, you can link the performance of your adherent cell line with the advantages of stirred-tank bioreactors. Fibra-Cel is a solid support growth matrix for mammalian and insect cells, used predominantly for production of secreted products such as recombinant proteins and viruses. When combined with the Eppendorf proprietary packed-bed impeller system, Fibra-Cel Disks enable growth of attachment-dependent cells in stirred-tank bioreactors and eliminate the need for cell filtration to separate cells from the end product.

## Performance


- > Proven culture success with a variety of cell lines
- > High surface-to-volume ratio supports high cell densities and productivities
- > Low shear environment facilitates the cultivation of sensitive cells

## Flexibility

- > Usable in glass, single-use, and sterilize-in-place bioreactors
- > Large working volume range of compatible bioreactors: 1.9 L to 32 L

## Ease of Use

- > Single-step inoculation
- > Easy product harvest
- > Easy setup of perfusion processes





# Performance


Powerful tools for your success

By providing a growth support, you can make the cultivation of anchorage-dependent cells in stirred-tank bioreactors a reality. Fibra-Cel Disks consist of a three-dimensional mesh to which cells attach, which protects them from damaging shear forces, and which simplifies the harvest of secreted products. Like this, Fibra-Cel Disks can help you in establishing highly efficient cell culture bioprocesses; in terms of cell and product yield as well as workload.




## Made to be populated

- > Attachment matrix with a diameter of 6 mm
- > Electrostatically pretreated to support cell adhesion



## Micro-structured

- > Made of a polyester/polypropylene mesh
- > High surface-to-volume ratio



## Growth-supporting

- > Provide a low-shear environment, as the cells become embedded in the three-dimensional fiber system

## Application examples\*

### A selection of cell types successfully grown on Fibra-Cel Disks

- > Vero ([1](#), [2](#), [11,13](#))
- > Adherent HEK 293 ([3](#), [4](#), [5](#), [6](#))
- > BHK ([11](#))
- > CHO ([7](#))
- > Mesenchymal stem cells ([8](#))
- > Hybridoma ([9](#))
- > Insect cells ([12](#))

### Typical products and applications

- > Recombinant protein expression ([4](#), [5](#), [6](#), [7](#), [9,12](#))
- > rVSV production ([1](#))
- > Lentiviral vector production ([3](#))
- > Rabies virus production ([10](#), [11](#))
- > Stem cell expansion ([8](#))

\*Please find the reference list on the last page of this brochure. The examples and references are non-comprehensive.



# Performance

Powerful tools for your success

By combining Fibra-Cel Disks with the Eppendorf proprietary packed-bed impeller system, you create a three-dimensional matrix. A low pressure drop across the bed minimizes variability and maintains viability of cells over the entire bed. The packed-bed technology convinces with a high and uniform mass transfer of nutrients and oxygen.

## The packed-bed impeller is build into the bioreactor to form a packed-bed bioreactor:

The packed-bed impeller forms a compartment inside the bioreactor, in which the Fibra-Cel Disks and the cells attached to them are entrapped.

### Impeller shaft:

Rotating shaft creates an under-pressure, which causes the medium to flow through the shaft and is expelled at the top.



### Fibra-Cel Disks:

Growth matrix filled into the packed-bed impeller



### Bioreactor headplate:

Can be equipped with sensors and tubes for liquid addition and withdrawal

### Packed-bed:

Culture medium and secreted products can pass freely in and out the bed

### Gassing:

Aeration of the medium takes place inside the impeller tube, providing a bubble-free aeration in and around the packed-bed, which is advantageous for the cultivation of shear sensitive cells.



# Flexibility

Use the bioreactor size and material of your needs

You can use Fibra-Cel Disks in combination with glass, single-use, and sterilize-in-place bioreactors of different sizes; dependent on your needs. With Eppendorf packed-bed bioreactors you can realize a growth area of up to  $1.4 \times 10^6 \text{ cm}^2$ , which offers you an enormous potential for process scale-up.

## Glass vessels

The Eppendorf glass vessels equipped with packed-bed impeller offer you great flexibility in working volumes; from a bed volume of 0.5 L to 5 L, corresponding to a growth surface of  $6 \times 10^4 \text{ cm}^2$  to  $6 \times 10^5 \text{ cm}^2$ .



### BioBLU® 5p Single-Use Bioreactor

The BioBLU 5p Single-Use Bioreactor combines the advantages of single-use vessels with the performance of packed-bed technology. The BioBLU 5p is preloaded with Fibra-Cel Disks - ready to use, saving you time on setup and cleaning.

The BioBLU 5p Single-Use Bioreactor has a bed volume of 1.5 L, corresponding to a growth surface of  $1.8 \times 10^5 \text{ cm}^2$ .



## Sterilize-in-place bioreactors

For higher capacity, a sterilize-in-place vessel equipped with a packed-bed impeller is available. The bed volume of is 12.0 L, corresponding to a surface area of  $1.4 \times 10^6 \text{ cm}^2$ .



## Overview of available vessel sizes

Working volume	Bed volume	Vessel type	Fibra-Cel Disks per vessel <sup>1</sup>	Total disk surface area	Roller bottle equivalent <sup>2</sup>
1.9 L	0.5 L	Glass	50 g	$6.0 \times 10^4 \text{ cm}^2$	71
3.8 L	1.5 L	Glass and single-use	150 g	$1.8 \times 10^5 \text{ cm}^2$	212
5.6 L	2.5 L	Glass	250 g	$3.0 \times 10^5 \text{ cm}^2$	353
10.5 L	5.0 L	Glass	500 g	$6.0 \times 10^5 \text{ cm}^2$	706
32.0 L	12.0 L	Stainless steel	1,200 g	$1.4 \times 10^6 \text{ cm}^2$	1,694

<sup>1</sup> On average, 10 g Fibra-Cel Disks are used per 100 mL bed volume

<sup>2</sup> The roller bottle equivalent is based on an average surface area of 850  $\text{cm}^2$  per bottle

eppendorf

Constant P/V-Based Scale-Up of Packed-Bed Perfusion Bioreactors for Culture of Adherent Cells



## Bioprocess scale-up in packed-bed bioreactors

Download this poster to discover, how to determine packed-bed vessel power numbers of bench-scale single-use bioreactors and production-scale sterilize-in-place (SIP) bioreactors to enable constant power per volume (P/V) scale-up.

[www.eppendorf.group/poster-fibra-cel](http://www.eppendorf.group/poster-fibra-cel)





# Ease of Use


Simple procedures from inoculation to harvest

When you inoculate your bioreactor, the cells become entrapped within the Fibra-Cel Disk bed, where they remain throughout the process. Fibra-Cel Disks enable sustained long-term periods of high-density growth in perfusion, without danger of clogging. They also simplify downstream processing, as they save you the need for subsequent cell filtration to separate cells from the end product.

## Efficient cell attachment

Cells usually attach efficiently to the Fibra-Cel Disks, which allows you inoculation in a straightforward manner.

- > Cell attachment takes usually only up to one hour
- > Cell attachment is efficient, because cells become imbedded in the fiber system and cannot fall off
- > No need for periodically stopping stirring to allow time for cells to become attached



## Easy product harvest

As the cells adhere to the Fibra-Cel Disks, they remain in the packed-bed throughout the bioprocess. Therefore you can easily harvest medium and secreted products, without the need for cell separation.

## Easy setup of perfusion processes

The packed-bed technology allows you to easily setup a perfusion processes. Fresh medium can be added to the bioreactor through a headplate port. Used medium can be removed from the bioreactor while the cells remain in the packed-bed. This makes filtration unnecessary and therefore avoids the danger of filter clogging.



## Vero cell culture in perfusion

Watch our webinar to discover, how a Vero cell perfusion process was set up using a packed-bed bioreactor and how it facilitated achieving high cell densities.

[www.eppendorf.group/webinar-fibra-cel](http://www.eppendorf.group/webinar-fibra-cel)



## Ordering Information

### Ordering information

Description	Order no.
<b>Fibra-Cel® Disks, 50 g</b>	M1292-9984
<b>Fibra-Cel® Disks, 150 g</b>	M1292-9992
<b>Fibra-Cel® Disks, 250 g</b>	M1292-9988
<b>Fibra-Cel® Disks, 1 kg</b>	M1292-9974

### Ordering information

Description	Order no.
<b>BioBLU® 5p Single-Use Bioreactor</b> , cell culture, microsparger, packed-bed impeller, optical pH, sterile, 1 piece	M1363-0119
<b>BioBLU® 5p Single-Use Bioreactor</b> , cell culture, macrosparger, packed-bed impeller, optical pH, sterile, 1 piece	M1363-0133

For more information please visit [www.eppendorf.com/Fibra-Cel](http://www.eppendorf.com/Fibra-Cel) and [www.eppendorf.com/BioBLUc](http://www.eppendorf.com/BioBLUc)

## Literature

- [1] Makovitzki et al. Evaluation of a downstream process for the recovery and concentration of a cell-culture-derived rVSV-Spike COVID-19 vaccine candidate. *Vaccine*, 39(48): 7044-7051, 2021
- [2] Guertin. A packed-bed bioreactor system for enhancing Vero cell growth in a semi-continuous mode of operation. *Bachelor thesis*, Worcester Polytechnic Institute, 2016
- [3] McCarron et al. Transient lentiviral vector production using a packed-bed bioreactor system. *Human Gene Therapy Methods*, 30(3), 2019
- [4] Dalton et al. Over-expression of secreted proteins from mammalian cell lines. *Protein Science*, 23(5): 517-525, 2014
- [5] Kaufman et al. Continuous production and recovery of recombinant Ca<sup>2+</sup> binding receptor from HEH 293 cells using perfusion through a packed bed bioreactor. *Cytotechnology*, 33: 3-11, 2000
- [6] Kim et al. MKR mice are resistant to the metabolic actions of both insulin and adiponectin: discordance between insulin resistance and adiponectin responsiveness. *American Journal of Physiology-Endocrinology and Metabolism*, 291(2): E298-E305, 2006
- [7] Hatton. Productivity studies utilizing recombinant CHO cells in stirred-tank bioreactors: A comparative study between the pitch-blade and packed-bed bioreactor systems. *Master thesis*, Utha State University, 2012
- [8] Tsai et al. Expansion of human mesenchymal stem cells in fibrous bed bioreactor. *Biochemical Engineering Journal*, 108: 51-57, 2016
- [9] Golmakany et al. Continuous production of monoclonal antibody in a packed-bed bioreactor. *Biotechnology and Applied Biochemistry*, 41(3): 273-278, 2010
- [10] Hassanzadeh et al. High Vero cell density and Rabies virus proliferation on Fibra Cel Disks versus Cytodex-1 in spinner flasks. *Pakistan Journal of Biological Sciences*, 14(7): 441-448, 2011
- [11] Gümüşderelioğlu et al. Rabies virus production in non-woven polyester fabric packed-bed bioreactors. *Biotechnology and Applied Biochemistry*, 33(3): 167-172, 2001
- [12] Kompier et al. Use of a stationary bed reactor and serum-free medium for the production of recombinant proteins in insect cells. *Enzyme and Microbial Technology*, 13(10): 822-827, 1991
- [13] Han et al. High-density Vero cell perfusion culture in BioBLU® 5p Single-Use Vessels. *Eppendorf Application Note 359*, 2017

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com/bioprocess](http://www.eppendorf.com/bioprocess)