

Eppendorf Certificate

Certificate of Quality – epT.I.P.S.[®] Motion Quality Assurance

All epT.I.P.S.[®] Motion pipette tips are single-use plastic tips. In combination with all Eppendorf epMotion[®] automated liquid handling systems, they form a perfectly functioning pipetting system.

Quality Assurance

The Eppendorf Quality Management System is certified according to the international standards ISO 9001:2015-09, ISO 14001:2015-09, and ISO 13485:2016-03.

All epT.I.P.S.[®] Motion pipette tips are produced in a controlled environment according to ISO class 8 of ISO 14644-1 and are subject to regular quality inspections within a statistical process control.

Composition

epT.I.P.S.[®] Motion pipette tips are made of virgin polypropylene of highest purity and quality. Our material suppliers do not use or intentionally incorporate the following agents into the materials Eppendorf uses for production:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyl dodecyl ammonium salts (DiHEMA))
- Plasticizers (softeners)
- Latex

Eppendorf confirms that all plastic materials fulfill the requirements as described in the current version of the FDA, CFR, Title 21 (Food and Drugs) in chapter 178.2010 “Antioxidants and Stabilizers for Polymers” and chapter 177.1520 (a)(1), (b) and (c)1.1 “Olefin polymers” regarding substances used for manufacturing of materials and articles or components of articles intended to come into contact with food.

epT.I.P.S.[®] Motion pipette tips and all product components do not contain:


- Metallic dyes
- Glucan
- Cellulose

All resin material is tested, qualified, and shows non-toxic biocompatibility according to the Standards USP Class VI Chapter 88, USP 661.1., EP 3.1.3, EP 3.1.6, and 10993-4:2017-12, ISO 10993-5:2009-10, ISO 10993-10:2020-04, ISO 10993-11:2018-09.

Hamburg, April 2024



Joana Tziolis
Product Life Cycle Manager
Division Consumables



Thomas Hengstmann
Head of Global Quality Operations

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.[®], are registered trademarks of Eppendorf SE, Germany.
All rights reserved, incl. graphics and images. Copyright ©2024 by Eppendorf SE.